# SLEslide

## ANTI-NATIVE DEOXYRIBONUCLEOPROTEIN ANTIBODIES

### Intended use:

Rapid Slide test for the qualitative and quantitative determination of anti-native deoxyribonucleoprotein (anti-DNP) antibodies in human serum and plasma.

### Summary:

The presence of autoantibodies to nuclear proteins is a common finding in systemic erythematosus (SLE) and other collagen diseases. Test currently used in the clinical laboratory for the detection of various antinuclear antibodies include the SLE tests, latex agglutination, complement fixation and immunofluorescence. Of the antinuclear antibodies, anti – desoxyribonucleoprotein (anti-DNP) appears to be the most specific for SLE.

Anti-DNP is present in high titers in the serum of the majority of SLE patients with active disease but is present only occasionally in remission states. Although anti-DNP is not found exclusively in SLE, only low titers may be detected in diseases such as chronic hepatitis, periarteritis nodosa, dermatomyositis, scleroderma and drug hypersensitivity. Biocon SLEslide is based on the principle of the latex agglutination assay described by Singer.

### Principle:

Biocon SLEslide is based on an immunologic reaction between DNP bound to biologically inert latex particles and anti-DNP in the test specimen. When serum containing anti-DNP is mixed with the latex reagent, visible agglutination occurs.

### Reagent concentration:

#### R1:

Polystyrene Latex particles coated with DNP prepared from calf thymus

Preservative

#### R2/positive control:

Antigen containing solution

Preservative

#### R3/negative control:

Solution without Antigen

Preservative

### Preparation and Stability:

The reagent and control sera are stable up to the expiry date when stored at +2 to +8°C. Do not freeze!

The reagent must be resuspended gently but completely prior use. We recommend use of a rolling mixer for 2 minutes.

Avoid air bubbles.

### Specimen:

Use fresh clear serum or plasma specimens.

Lipemic serum or microbial contamination may cause erroneous results. If the test cannot be performed immediately, store the specimen at +2 to +8°C for up to 24 hours. For longer storage, freeze the serum.

### Notes:

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Reagents containing sodium azide may combine with copper and lead plumbing to form highly explosive metal azides.

Dispose of reagent by flushing with large amounts of water to prevent azide buildup.

The positive and negative controls were prepared from human sera which have been tested using standard methods and found to be non-reactive for HbsAg and HIV antibodies. However, no test method can offer complete assurance that infectious agents are absent. Therefore all human specimens should be considered potentially infectious.

### Limitations - interference:

Anti-DNP is usually present in active SLE and occasionally present in inactive SLE. A number of drugs cause an SLE like syndrome accompanied with positive anti-DNP titer. The symptoms and titer disappear upon drug withdrawal

For diagnostic purposes, SLEslide results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Anti-DNP may be found in diseases other than SLE. Low titers have been detected in rheumatoid arthritis, chronic hepatitis, periarteritis nodosa, dermatomyositis, scleroderma, atypical, pneumonia, tuberculosis and lymphoma.

Reaction time is critical. If the reaction exceeds 3 minutes, drying of the reaction mixture may cause false positive results.

Freezing the SLE Latex Reagent will result in spontaneous agglutination.

Intensity of the agglutination is not necessarily indicative of relative SLE titers; therefore screening reactions should not be graded.

### Testing procedure:


#### Qualitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Place one drop (1 x 50 µl) of the positive control on field # 1 of the reaction slide.
3. Place one drop (1 x 50 µl) of the negative control on field # 2. Using a pipette, place one drop (50 µl) of each test specimen on successive fields.
4. Gently resuspend the Latex Reagent and add one drop (35 µl) to each test field. Use stick to spread reaction mixture over entire test field.
5. Rotate the slide for 3 minutes and read **immediately** under direct light.

#### Quantitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Using saline, dilute the specimens 1:2, 1:4, 1:8, 1:16 or as needed.
3. Place one drop (50 µl) of each dilution on successive fields of the reaction slide.
4. Resuspend the SLEslide Latex Reagent gently and add one drop (50 µl) to each test field. Use stick to spread reaction mixture over entire test field.
5. Rotate the slide for 3 minutes and read **immediately** under direct light.

### Interpretation:

### Results:

#### Qualitative test:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the negative control (2). A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared with the positive control (1).

#### Quantitative test:

The titer of the serum is the reciprocal of the highest dilution which exhibits a positive reaction.

# SLEslide

## ANTI-NATIVE DEOXYRIBONUCLEOPROTEIN ANTIBODIES

### Sensitivity and method comparison:

Serum specimens from 210 individuals were examined. The Biocon SLEslide, another commercial SLEslide test and in 145 cases an additional immunofluorescence test were performed. In the case of ANA titers of 1:10 no signal was detected. Test results with higher titers show no significant difference between both Latex assays. These data indicate that both sensitivity and specificity of the Biocon SLEslide Latex are 100%.

### Precision:

In a study on precision, a panel of 10 serum samples with SLE titers ranging from 1.1 to 1.64 were tested on 10 consecutive days by the Quantitative Method (100 determinations). No determinations gave more than a 2-fold difference from the mean titer of a sample.

### Quality control:

Biocon SLEslide **R2/positive control** and **R3/negative control** should be included in each test series. The **R2** should produce strong agglutination and the **R3** should yield a smooth suspension with no agglutination.

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2. Hughes, G.R.V., Cohen, S.A. and Christian, C.L., Ann. Rheum. Dis., 30, 259-264 (1971).